



Radiation Risk Assessment



READ THE NOTES ON COMPLETION BEFORE ATTEMPTING TO FILL IN THIS FORM. IF INSUFFICIENT SPACE IS AVAILABLE UNDER ANY SECTION, USE A SEPARATE PIECE OF PAPER AND ATTACH IT TO THE FORM.			FILE REF: DATE: 6 th May 2011		
DEPARTMENT: Cardiovascular Medicine Wellcome Trust Centre for Human Genetics		PERSONS INVOLVED: Bhattacharya Group See attached sheet		OTHERS AT RISK: None	
LOCATION OF WORK: Room 10/088 (Lab 3)					
DESCRIPTION OF PROCEDURE: <h2 style="text-align: center;">Hybridisation of Southern Blots with ³²P – see attached protocols</h2>					
SUBSTANCES USED	QUANTITIES USED	FREQUENCY OF USE	HAZARDS IDENTIFIED	EXPOSURE ROUTE	DOSE PER PROCEDURE (µSv)
[α-P32]dCTP	1.85 MBq or 50µCi per reaction (taken from stock pot containing max 9.25 MBq)	Less than once per month	High energy β-emitter (1.71MeV). Skin dose rate from 1MBq point source at 30cm is 0.12mSv/h [Delacroix et al. <i>Radiat. Prot. Dosim.</i> 98 , 2002]	Skin (mostly fingers)	<i>Whole body:</i> Effectively zero <i>Extremities:</i> 300µSv per probe
COULD A LESS HAZARDOUS SUBSTANCE (OR FORM OF THE SUBSTANCE) BE USED INSTEAD? Yes/No JUSTIFY NOT USING IT:					
COULD A LOWER ACTIVITY BE USED? Yes / No JUSTIFY QUANTITY OF MATERIAL IN USE: Quantities used are specified in the probe labelling kit.					
WHAT MEASURES HAVE YOU TAKEN TO CONTROL RISK? ENGINEERING CONTROLS: Perspex screen and radiation protection cabinet. Stock pot shipped and stored encased in lead. Tubes kept in Perspex racks with lids. Hybridisation tubes containing labelled DNA carried from 10/088 to hybridisation oven in lidded Perspex box. Reduce dosage by control of distance, time and shielding.					
PPE: Lab coat Nitrile gloves Safety goggles					
MANAGEMENT MEASURES: Follow the As Low As Reasonably Practicable (ALARP) rule Adhere to local rules, EPR2010 and IRR99 Monitor work area before, during and after use Lone working prohibited Work only in designated Supervised areas Radioactivity stock pot stored in locked refrigerator in Room 10/088 Adhere to limits of designated sink and bins Adhere to storage time limits for bins [Green/yellow band bin - 6 months; Orange bin - 6 months]					
CHECKS ON CONTROL MEASURES: Regular monitoring and supervision Checks on documentation – control of Monthly returns					
Radiation monitor: EP15		TRAINING REQUIREMENTS: URPO lecture “Working with unsealed radioactive sources” and in-house training.			
Is dosimetry required? Yes / No Type of dosimetry: Extremity / Whole body / Biological					
EMERGENCY PROCEDURES: Refer to Contingency Plans in Local Rules & University Policy Statement S8/05: Appendix 16. Decontamination where necessary with Count-Off and/or Decon90 decontaminants as recommended. Dispose of contaminated materials in designated bins and sinks. For any incident beyond minor contamination of radioactive workspaces alert RPS (and SRPS/URPO if appropriate).			WASTE DISPOSAL: Aqueous ~50% to designated sink Solid ~50% to green-yellow band bin Actual % amounts estimated after DNA labelling using EP15 or similar monitor		

Are overall risk control measures adequate? Yes/No	SIGNATURE
NAME AND POSITION OF ASSESSOR:	SIGNATURE:
NAME OF RADIATION PROTECTION SUPERVISOR: Dr James Brown	SIGNATURE:

Date of routine review	DATE:	/ / 2012	/ / 2013	/ / 2014	/ / 2015	/ / 2016
	BY:					



Department of Cardiovascular Medicine
University of Oxford
Radiation Risk Assessment Training Record



Title: Hybridisation of Southern Blots with ^{32}P

The undersigned have read the above radiation risk assessment and understand the safety arrangements required and their own obligations in ensuring compliance.

DATE	NAME & SIGNATURE	TRAINED BY (NAME & SIGNATURE)

Phosphorus - 32

³²P₁₅






Half life: 14.3 days
Specific activity: 1.06E+16 Bq.g⁻¹

Risk group: 2
Risk colour: Orange



Main emissions (keV)					Exemption levels			
	Gamma or X		Beta (E _{max})		Electrons		Alpha	
	E	%	E	%	E	%	E	%
E1			1710	100				
E2								
E3								
% omitted			0					

Exemption levels	
Quantity (Bq)	1E+05
Concentration (Bq.g ⁻¹)	1E+03

Transport (TBq)	
IAEA ST1 A ₁ value	5E-1
IAEA ST1 A ₂ value	5E-1

EXTERNAL EXPOSURE (mSv.h ⁻¹) for an activity of 1 MBq or 1 MBq.m ⁻² (as appropriate)				
Point source (30 cm)	Infinite plane source	10 ml glass vial	Contact with 50 ml glass beaker	Contact with 5 ml plastic syringe
				
<i>Betas, electrons (skin dose)</i>	<i>Betas, electrons (skin)</i>			
1.18E-1	10 cm: 1.4E-01 1 m: 4.8E-02			
<i>Gammas, X rays (deep tissue dose)</i>	<i>Photons (skin)</i>			
0.00E+0	10 cm: 0.0E+00 1 m: 0.0E+00	100 cm: 1.31E-6	7.11E-4	2.39E+1
	<i>Photons (deep dose)</i>			
	10 cm: 0.0E+00 1 m: 0.0E+00			

The values above do not include Bremsstrahlung radiation.

CONTAMINATION			SHIELDING (mm)																			
Contamination skin dose (mSv.h⁻¹)	Detection	Derived limits (Bq.cm⁻²)	Betas and electrons (Total absorption)																			
Uniform deposit (1kBq.cm ⁻²): 1.89E+0	<table border="1"> <tr><th colspan="2">Recommended probes*</th></tr> <tr><td>Alpha</td><td></td></tr> <tr><td>Beta</td><td>++</td></tr> <tr><td>Gamma</td><td></td></tr> <tr><td>X rays</td><td></td></tr> </table>	Recommended probes*		Alpha		Beta	++	Gamma		X rays		<table border="1"> <tr><th colspan="2">Removable contamination</th></tr> <tr><td>5E+1</td><td></td></tr> <tr><th colspan="2">Fixed contamination</th></tr> <tr><td>3E+2</td><td></td></tr> </table>	Removable contamination		5E+1		Fixed contamination		3E+2		Glass: 3.4	Plastic: 6.3
Recommended probes*																						
Alpha																						
Beta		++																				
Gamma																						
X rays																						
Removable contamination																						
5E+1																						
Fixed contamination																						
3E+2																						
0.05 ml droplet (1 kBq): 1.33E+0			Gamma and X rays (half and tenth value thickness)																			
				1/2	1/10																	
			Lead: -	Steel: -																		

* If no probes are indicated the recommended technique is to use a wipe test in association with a probe or liquid scintillation technique

INTERNAL EXPOSURE FOR WORKERS					
COMMITTED EFFECTIVE DOSE PER UNIT INTAKE (Sv.Bq ⁻¹)					
Ingestion	f₁		Inhalation	1 μm	5 μm
All compounds	0.800	2.4E-09	All unspec. compounds	F: 8.0E-10	1.1E-09
			Some phosphat. det. by comb. cation	M: 3.2E-09	2.9E-09
				S	
Highest dose organ	Lungs	20 mSv A _{LI,ingestion}	8.3E+06 (Bq)	20 mSv A _{LI,inhalation}	6.3E+06 (Bq)

MAXIMUM RECOMMENDED ACTIVITIES IN LOW LEVEL OR INTERMEDIATE LEVEL LABORATORIES (Bq)							
PHYSICOCHEMICAL STATE	Subject to external exposure requirements which may be more restrictive						
	Volatility factor (k)	Supervised area			Controlled area		
		Bench	Fume hood		Bench	Fume hood	Glove box
All compounds	0.01	5E+05	5E+06	2E+06	2E+07	2E+09	

Dose estimation - Hybridisation of Southern Blots with ³²P

ALMOST ALL WORK WILL BE CARRIED OUT EITHER BEHIND A PERSPEX SHIELD OR USING PERSPEX TUBE RACKS SO THE DOSE SHOULD BE LIMITED TO FINGERS; WHOLE BODY DOSE WILL BE INSIGNIFICANT.

<i>Handling stock pot</i>	Maximum activity 9.25MBq
<i>Radiolabel used per probe</i>	Maximum activity of 1.85MBq
<i>Assume 50% incorporation (varies between probes)</i>	Maximum activity of 0.925MBq per labelled probe
<i>Washing away unhybridised probe (progressively lower activity handled)</i>	Maximum activity of 0.0148MBq per paper

Dose estimation per probe (some experiments may use up to 3 probes)

	For an activity of 1MBq, from Delacroix				dose (mSv)		
	High dose	Syringe	Low dose	Point source (30cm)			
	H	23.9 mSv/hr	L	0.12 mSv/hr			
		0.398333333 mSv/min		0.002 mSv/min			
		0.006638889 mSv/sec		3.33333E-05 mSv/sec			
	<i>dose rate</i>	<i>manipulation</i>	<i>manipulated</i>				
	<i>enter H or L</i>	<i>time (s)</i>	<i>amt (MBq)</i>				
<i>Pipetting from stock pot into probe labelling reaction</i>	H	3	9.25	0.1842			
<i>Pipetting up and down</i>	L	10	1.85	0.0006			
<i>Closing lids and transfer to/from centrifuge</i>	H	5	1.85	0.0614			
<i>Adding probe to Quick Spin column for clean-up</i>	L	5	1.85	0.0003			
<i>Centrifugation & disposal of column</i>	H	5	1.85	0.0614			
<i>Adding to hybridisation buffer (assuming 50% incorporation)</i>	L	10	0.925	0.0003			
<i>Pouring off hybridisation buffer</i>	L	10	0.925	0.0003			
<i>Wash 1</i>	L	20	0.185	0.0001			
<i>Wash 2</i>	L	20	0.0185	0.0000			
<i>Wash 3 etc</i>	L	20	0.00185	1.23333E-06			
			Total	0.3087	Total finger dose in mSev		

Southern Blotting Protocol

Prepare the gel

Digest DNA samples with appropriate restriction enzyme(s).

The amount of DNA per well for hemizygous transgenic mice should be 10µg. Prepare copy number standards of the transgene to construct a standard curve to determine transgene copy number from unknown samples.

Make an agarose gel. Use 1X TAE, and pour the gel no more than 0.7cm thick. Add loading dye to digests.

Load digested samples onto 0.8% TAE gel without EtBr (*EtBr can alter the mobility of DNA, confounding measurement, so stain afterwards*):

Electrophorese relatively slowly until the dye passes off the end of the gel. (*Running the gel slowly results in thinner, sharper bands on the Southern*).

EtBr stain and photograph gel with ruler from bottom of wells

Wash the Gel

Rinse stained gel briefly in dH₂O

DEPURINATION – Wash in 1L 0.12M HCl for 10 min on rocker, then rinse in dH₂O

DENATURATION – 30 mins in 1L 1.5M NaCl/ 0.5M NaOH, then rinse in dH₂O

Set up the transfer

Alkalyne transfer: DNA is transferred to Hybond XL membranes by capillary blotting in denaturation buffer.

Set up: Cling wrap, gel (turned over), Hybond XL membrane (parafilm the edges to prevent “short circuit”), 3-5 sheets of Whatman paper, approximately 2 inches of paper towels, flat plastic and weight down.

Transfer O/N

Disassemble the transfer pyramid

Recover the membrane & mark position of wells in pencil and ensure that the up-down and back-front orientations are recognisable.

Immobilise the DNA

Irradiate in UV crosslinker at 70,000 µJ/cm²

Store membrane between sheets of 3MM for later hybridisation

Use blunt forceps to handle membrane

PROTOCOL INVOLVING THE USE OF RADIOACTIVE ISOTOPES

- PROCEDURE:** Hybridisation of Southern blots with ^{32}P
- ISOTOPE:** $\text{P}^{32}\alpha\text{-dCTP}$ (Redivue)
- STORAGE:** Locked Box in $+4^\circ\text{C}$ refrigerator in room 10/088
- QUANTITY:** 1.85 MBq per reaction

Prehybridisation

Pre-warm *Modified Church & Gilbert Buffer* (Prehybridisation Buffer):

7% (w/v) SDS

0.5 M Phosphate buffer, pH 7.2 (Na_2HPO_4)

10 mM EDTA

1% BSA

Pre-wet membrane in 6xSSC

Prehybridise membrane in 15mL hybridisation buffer (DNA side facing in) for 1 hr at 65°C

Prepare Probe

Dilute 25ng probe in 45 μL of 10/1 Tris/EDTA pH8.0

Denature the DNA at $95\text{-}100^\circ\text{C}$ for 5 mins

Snap cool on ice for 5 minutes

Centrifuge briefly at max speed

Labelling Probe

Label using Amersham Rediprime II Random Prime Labelling system (GE Healthcare, RPN1633)

(In hot suite)

Add denatured probe to the reaction tube (do NOT mix at this stage)

Add 5 μL $\text{P}^{32}\alpha\text{-dCTP}$ to denatured probe and mix by pipetting up and down about 12 times, moving the pipette tip around in the solution.

Incubate at room temperature for at least 30 minutes.

Purification of probe using Sephadex G50 spin columns

Use Roche mini Quick Spin Columns #11 814 419 001

Preparing the column

1. Resuspend column matrix by flicking column sharply several times
2. Remove top cap, then snap off bottom tip, and place column into sterile microfuge tube.
3. Spin column (1000g x 1 min) to pack the column and remove residual buffer

Purifying the sample

Place prepared column in clean, sterile microfuge tube (preferably screw top)

Carefully apply sample to centre of column bed

Spin column (1000 g, 4 min)

Quickly check % incorporation by comparing counts in eluate vs column

Discard column

Hybridisation of labelled probe to the blot

Denature probe again at 95-100°C for 5 mins, then snap cool on ice for 5 mins

Change the hybridization buffer in the tube. Allow to warm up to 65 degrees.

Add denatured probe directly to the hybridisation buffer already in the hybridisation tube.

Hybridise overnight at 65°C

Next day. Washing the blot

Pour off liquid from hybridisation bottle straight down the designated sink and flush sink with tap water.

Carry out a series of washes, of differing stringency depending on the probe:

Rinse briefly in 2X SSC/0.1% SDS

Standard wash steps:

2X SSC/1% SDS, 65 degrees, 30 mins

0.5X SSC/0.1% SDS, 65 degrees, 15 mins

0.2X SSC/0.1% SDS, 65 degrees, 10 mins

Remove blots from the hybridisation tubes and rinse briefly in either low stringency buffer or 2XSSC

Wrap blot in Saran wrap and autoradiograph.

WASTE:

Depending on % incorporation

Liquid disposed activity: % in probe (~50%)

Solid disposed activity: % remaining in Sephadex column (~50)

- Record all solid and liquid waste on the appropriate disposal forms.
- Monitor work area, disposal sink and hybridisation oven for contamination.
- Fill in the relevant monitoring sheets.