




Radiation Risk Assessment



				DATE: 21 st Aug 2017	
DEPARTMENT: Cardiovascular Medicine Wellcome Trust Centre for Human Genetics			PERSONS INVOLVED: See attached training record		OTHERS AT RISK: None
LOCATION OF WORK: Room 10/088 (Lab 3)					
DESCRIPTION OF PROCEDURE: ¹⁴C-creatine uptake by cells – see attached protocols.					
SUBSTANCES USED	QUANTITIES USED	FREQUENCY OF USE	HAZARDS IDENTIFIED	EXPOSURE ROUTE	DOSE PER PROCEDURE (µSv)
¹⁴ C-creatine hydrate	Between 37 and 74 kBq (1-2 µCi) per experiment	2-3 times per week	Low energy β-emitter (157 keV). Skin dose rate from 1MBq point source at 30cm is 0 mSv/h. Bremsstrahlung radiation may be significant. [Delacroix et al. <i>Radiat. Prot. Dosim.</i> 98 , 2002]	Skin (mostly fingers)	Whole body: Effectively zero
					Extremities: Effectively zero
COULD A LESS HAZARDOUS SUBSTANCE (OR FORM OF THE SUBSTANCE) BE USED INSTEAD?				Yes/No	
COULD A LOWER ACTIVITY BE USED?			No		
JUSTIFY QUANTITY OF MATERIAL IN USE: Quantities used are already very low and cannot be further reduced due to sensitivity levels required in assay.					
WHAT MEASURES HAVE YOU TAKEN TO CONTROL RISK?					
ENGINEERING CONTROLS: Work carried out in designated areas on radioactivity spill trays. Reduce dosage by control of distance, time and shielding.					
PPE: Lab coat, nitrile gloves and safety eyewear					
MANAGEMENT MEASURES: Follow the As Low As Reasonably Practicable (ALARP) rule Adhere to local rules, EPR2010 and IRR99 Monitor work area before, during and after use. Wipe tests taken after each set of experiments. Lone working prohibited Work only in designated Supervised areas Radioactivity stock pot stored in locked refrigerator in Room 10/088 Adhere to limits of designated sink and bins Adhere to storage time limits for bins [12 months for both 'solids' green/yellow band bin & 'organic liquids' orange bin]					
CHECKS ON CONTROL MEASURES: Regular monitoring (including wipe tests by users and RPSs) and supervision Checks on documentation – control of Monthly returns					
Radiation monitor: EP15		TRAINING REQUIREMENTS: URPO lecture "Working with unsealed radioactive sources" and in-house training.			
Is dosimetry required? No					
EMERGENCY PROCEDURES: Refer to Contingency Plans in Local Rules & University Policy Statement S8/05: Appendix 16. Decontamination where necessary with Count-Off and/or Decon90 decontaminants as recommended. Dispose of contaminated materials in designated bins and sinks. For any incident beyond minor contamination of radioactive workspaces alert RPS (and SRPS/URPO if appropriate).			WASTE DISPOSAL: Approximations – actual amounts determined from scintillation counting Aqueous 89% to designated sink Solid 1% to green-yellow band bin Scintillant 10% to orange bin		

Are overall risk control measures adequate?	<p style="text-align: center;">Yes</p>
NAME OF RADIATION PROTECTION SUPERVISOR: <p style="text-align: center;">Dr James Brown</p>	SIGNATURE: 

Date of routine review	DATE:	/ / 2018	/ / 2019	/ / 2020	/ / 2021	/ / 2022
	BY:					

Carbon - 14



Half life: 5730 years
 Specific activity: 1.65E+11 Bq.g⁻¹

Risk group: 4
 Risk colour: Green

Main emissions (keV)								
	Gamma or X		Beta (E _{max})		Electrons		Alpha	
	E	%	E	%	E	%	E	%
E1			157	100				
E2								
E3								
% omitted			0					

Exemption levels	
Quantity (Bq)	1E+07
Concentration (Bq.g ⁻¹)	1E+04

Transport (TBq)	
IAEA ST1 A ₁ value	4E+1
IAEA ST1 A ₂ value	3E+0

EXTERNAL EXPOSURE (mSv.h ⁻¹) for an activity of 1 MBq or 1 MBq.m ⁻² (as appropriate)				
Point source (30 cm)	Infinite plane source	10 ml glass vial	Contact with 50 ml glass beaker	Contact with 5 ml plastic syringe
<i>Betas, electrons (skin dose)</i> 0.00E+0	<i>Betas, electrons (skin)</i> 10 cm 0.0E+00 1 m 0.0E+00	100 cm Brem. Rad.	Brem. Rad.	Brem. Rad.
<i>Gammas, X rays (deep tissue dose)</i> 0.00E+0	<i>Photons (skin)</i> 10 cm 0.0E+00 1 m 0.0E+00 <i>Photons (deep dose)</i> 10 cm 0.0E+00 1 m 0.0E+00			

The values above do not include Bremsstrahlung radiation. Brem. Rad. indicates that it may be significant.

CONTAMINATION												
Contamination skin dose (mSv.h⁻¹)	Detection	Derived limits (Bq.cm⁻²)										
Uniform deposit (1kBq.cm ⁻²) 3.24E-1	<table border="1"> <tr><th colspan="2">Recommended probes*</th></tr> <tr><td>Alpha</td><td></td></tr> <tr><td>Beta</td><td>++</td></tr> <tr><td>Gamma</td><td></td></tr> <tr><td>X rays</td><td></td></tr> </table>	Recommended probes*		Alpha		Beta	++	Gamma		X rays		Removable contamination
Recommended probes*												
Alpha												
Beta		++										
Gamma												
X rays												
0.05 ml droplet (1 kBq) 2.70E-3		4E+2										
		Fixed contamination										
		4E+4										

* If no probes are indicated the recommended technique is to use a wipe test in association with a probe or liquid scintillation technique

SHIELDING (mm)		
Betas and electrons (Total absorption)		
Glass	0.2	
Plastic	0.3	
Gamma and X rays (half and tenth value thickness)		
	1/2	1/10
Lead	-	-
Steel	-	-

INTERNAL EXPOSURE FOR WORKERS				
COMMITTED EFFECTIVE DOSE PER UNIT INTAKE (Sv.Bq ⁻¹)				
Ingestion	f₁	Inhalation (for soluble and reactive gases and vapours)		
Labelled organic compounds	1.000	5.8E-10	Vapour	5.8E-10
			Dioxide	6.2E-12
			Monoxide	8.0E-13
Highest dose organ	Whole body	20 mSv A _{LI} Ingestion	3.4E+07 (Bq)	20 mSv A _{LI} Inhalation
				3.4E+07 (Bq)

MAXIMUM RECOMMENDED ACTIVITIES IN LOW LEVEL OR INTERMEDIATE LEVEL LABORATORIES (Bq)							
PHYSICOCHEMICAL STATE	Subject to external exposure requirements which may be more restrictive						
	Volatility factor (k)	Supervised area			Controlled area		
		Bench	Fume hood		Bench	Fume hood	Glove box
All compounds except below	0.01	1E+07	1E+08	3E+07	3E+08	5E+09	
Oxide & anhydride	1	Forbidden	1E+06	Forbidden	3E+06	3E+08	

¹⁴C-creatine uptake by cells – Dose estimation

ALMOST ALL WORK WILL BE CARRIED OUT EITHER BEHIND A PERSPEX SHIELD OR USING PERSPEX TUBE RACKS SO THE DOSE SHOULD BE LIMITED TO FINGERS; WHOLE BODY DOSE WILL BE INSIGNIFICANT.

BREMSSTRAHLUNG RADIATION IS UNLIKELY TO BE SIGNIFICANT DUE TO THE LOW ENERGY OF THE RADIATION AND THE USE OF PERSPEX SHIELDING.


The external dose will be effectively zero since Perspex shielding, assay tubes, gloves and skin will all stop low energy beta-particles arising from ¹⁴C decay.

The maximum dose likely to be encountered arises from handling the stock material which is in screw-capped microcentrifuge tubes at a maximum activity of 4.6MBq (1250µl @ 3.7MBq/ml). This is handled for a matter of seconds and estimation of the dose is impractical based on the lack of dose values from Delacroix *et al.*

Radiolabelled Creatine Uptake Experiment Protocol

1. Plate cells into 6 well plates in preparation for protocol, including a duplicate plate which will not be spiked with isotope, but will be used for the determination of protein levels in each well.
2. Check hot lab for contamination by previous user prior to beginning work with the isotope. Can use Geiger counter, but using cotton buds dipped in **Decon 90** to wipe the area, then placing them in scintillation vials containing scintillation fluid is essential. These can be counted later together with everything else.
3. Prepare 13ml creatine-containing complete media (creatine concentration 500 μ M) in a 50ml Falcon tube. To this, add 10 μ l (37 kBq) 14 C-creatine.
4. Record the amount of isotope taken from the pot on the appropriate sheet.
5. Transport plate of cells to hot lab from the incubator, aspirate media and discard.
6. Add hot media containing creatine (2ml per well). Keep a volume of hot media (100 μ l) in order to measure later, to confirm the initial amount of isotope added to the cells.
7. Transport the 6 well plate containing the radiolabel-spiked media to the designated CO₂ incubator in 10/079 in a suitable Tupperware (or similar) container with a sealed lid.
8. Incubate the plate containing the radiolabelled creatine for 60 minutes alongside the duplicate (non 14 C spiked) plate of cells, which has standard complete media containing 500 μ M creatine, but no isotope. The hot plate should be placed on a suitable (ie lipped) tray in the incubator, in order that it would contain any spills if they were to occur. This plate should be clearly labelled as radioactive, and if possible, incubated on the bottom shelf of the incubator in order that spills cannot contaminate anything below it.
9. Return the plate to the hot suite in the sealed container. The duplicate plate can be dealt with on the bench and should not be taken to the hot suite.
10. Aspirate the hot media from the 6 wells. To do this, use disposable 3ml plastettes (disposable plastic Pasteur pipettes). Pipette the waste hot media into a 50ml Falcon tube.
11. Wash each of the 6 wells 3 times with 2ml PBS /well. This PBS and the media previously aspirated can also be placed in a 50ml Falcon tube, ready for disposal.
12. Solubilise cells with the addition of 0.5ml of 0.5% **Triton-X**.
13. Scrape cells free from the base of the wells using a cell scraper. This should be performed with the plate behind the Perspex shield, in order to eliminate the risk of any potential splash reaching you.
14. Aspirate the contents of each well using a P1000 micropipette with filter tips to a scintillation vial, then add a further 0.5ml of 0.5% Triton-X in PBS to the wells to rinse them, and finally aspirate this and add this also to the scintillation vials. Then, add 9ml **FluoranSafe XE** scintillation buffer to the vials, bringing the total volume in each one to 10ml. Also prepare 1 vial with 10ml of the scintillation buffer alone, and 1 with 9.9ml of the buffer and 100 μ l of the hot media which was saved at the beginning of the experiment.
15. Monitor the hot area using cotton buds or small pieces of filter paper dipped in Decon 90 to wipe the appropriate surfaces. Include pipette, Tupperware container etc as well as the area itself. Place buds in scintillation vials containing scintillation buffer, labelled with the area wiped.

16. Discard liquid waste collected in Falcon tubes down the designated sink, after first pouring a small volume of Decon 90 down the sink in order to aid the removal of the isotope. Record this on the sink waste sheet, with initials and details of isotope.
17. Discard the solid waste, such as tips and falcon tubes in the green solid waste bin, and mark the assumed amount of isotope on these items (eg 1% of total starting amount for tips).
18. Transport the vials to the scintillation counter in an appropriate sealed container, and count using carbon-14 program.
19. Discard scintillation waste in orange scintillation bin by the counter, and record it on the appropriate scintillation waste sheet (orange).
20. Finally, document the fate of all the waste on the original stock sheet where the amount of isotope taken was first recorded.



* - determined from scintillation counting